

# Certificate

## Thermo Scientific Nunc MaxiSorp Surface

This certificate is your guarantee that representative samples of this product have been evaluated and meet the following specifications:

**Test:**

Reproducibility/homogeneity of adsorption of IgG.

**Specification:**

Coefficient of variation (CV) within a plate, less than 5% well to well.  
OD results of individual wells are within +/- 10% of the plate mean.

**LOT**

Representative samples of this product produced by Nunc A/S, Thermo Fisher Scientific, have passed the tests in our Quality Control Department with satisfactory result.

In case of complaint, please fill in the form below and return it to your supplier.

Name: \_\_\_\_\_

Company: \_\_\_\_\_

Address: \_\_\_\_\_

Nunc distributor: \_\_\_\_\_

Date: \_\_\_\_\_ Signature: \_\_\_\_\_

Cause for the complaint:

Defective product enclosed:  Yes  No

Thermo Fisher Scientific · Nunc A/S  
Kamstrupvej 90 · P.O. Box 280  
DK-4000 Roskilde · Denmark  
T +45 4631 2000 · F +45 4631 2175  
E nunc.nalgene.eu@thermofisher.com

Thermo Fisher Scientific · Nalge Nunc International  
75 Panorama Creek Drive  
Rochester, NY 14625-2385 · U.S.A.  
T +1 585 586 8800 · F +1 585 586 8987  
E nunc.nalgene.na@thermofisher.com

[www.thermoscientific.com](http://www.thermoscientific.com)

© 2012 Thermo Fisher Scientific Inc. All Rights Reserved. - 71618 - Ver. 6.0 - 0512  
All trademarks are property of Thermo Fisher Scientific Inc. and its subsidiaries.

**Thermo**  
SCIENTIFIC  
Part of Thermo Fisher Scientific