

Certificate

Thermo Scientific™ Nunc™ MediSorp™ Surface

This certificate is your guarantee that representative samples of this product have been validated according to the following test.

Control of homogeneity in adsorption of IgG

Coefficient of variation (CV) within a plate format is less than 5%.
OD results in the wells of the plate are within $\pm 10\%$ from mean.

LOT

Representative samples of this product produced by Thermo Fisher Scientific, have passed the tests in our Quality Control Department with satisfactory result.

In case of complaint, please fill in the form below and return to your supplier.

Name: _____

Company: _____

Address: _____

Nunc distributor: _____

Date: _____ Signature: _____

Cause for the complaint:

Defective product enclosed: Yes No

Thermo Fisher Scientific · Nunc A/S
Kamstrupvej 90 · P.O. Box 280
DK-4000 Roskilde · Denmark
T +45 4631 2000 · F +45 4631 2175
E techsupport.nunc.eu@thermofisher.com

Thermo Fisher Scientific · Nalge Nunc International
75 Panorama Creek Drive
Rochester, NY 14625-2385 · U.S.A.
T +1 585 586 8800 · F +1 585 586 8987
E technicalsupport@thermofisher.com

www.thermoscientific.com

© 2013 Thermo Fisher Scientific Inc. All Rights Reserved. 71685 - Ver. 7.0 - 0913
All trademarks are property of Thermo Fisher Scientific Inc. and its subsidiaries.

Thermo
SCIENTIFIC
Part of Thermo Fisher Scientific