

Certificate

Thermo Scientific Nunc MaxiSorp Surface

This certificate is your guarantee that representative samples of this product have been evaluated and meet the following specifications:

Test:

Reproducibility/homogeneity of adsorption of IgG.


Specification:

Coefficient of variation (CV) within a lot, less than 5% tube to tube.
OD results of individual tubes are within +/- 10% of the lot mean.

This certificate concerns the following products:

Thermo Scientific Nunc Immuno Tube MaxiSorp

Application: Immunology



LOT

Representative samples of this product produced by Nunc A/S, Thermo Fisher Scientific, have passed the tests in our Quality Control Department with satisfactory result.

In case of complaint, please fill in the form below and return to your supplier.

Name: _____

Company: _____

Address: _____

Nunc distributor: _____

Date: _____ Signature: _____

Cause for the complaint:

Defective product enclosed: Yes No

Thermo Fisher Scientific · Nunc A/S
Kamstrupvej 90 · P.O. Box 280
DK-4000 Roskilde · Denmark
T +45 4631 2000 · F +45 4631 2175
E techsupport.nunc.eu@thermofisher.com

Thermo Fisher Scientific · Nalge Nunc International
75 Panorama Creek Drive
Rochester, NY 14625-2385 · U.S.A.
T +1 585 586 8800 · F +1 585 586 8987
E technicalsupport@thermofisher.com

www.thermoscientific.com

© 2013 Thermo Fisher Scientific Inc. All Rights Reserved.
All trademarks are property of Thermo Fisher Scientific Inc. and its subsidiaries.
74683 Ver. 6 - 0413

Thermo
SCIENTIFIC
Part of Thermo Fisher Scientific